# Orbegozo

TENSIÓMETRO DIGITAL DE MUÑECA – MANUAL DE INSTRUCCIONES  
WRIST BLOOD PRESSURE MONITOR- INSTRUCTION MANUAL  
TENSIO MÈTRE ELECTRONIQUE AU POIGNET - MANUEL D'INSTRUCTIONS  
TENSIÔMETRO DO PULSO DIGITAL- MANUAL DE INSTRUÇÕES

## TES 3650

Sonifer, S.A.  
Avenida de Santiago, 86  
30007 Murcia España  
[sonifer@sonifer.es](mailto:sonifer@sonifer.es)

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**European representative name:**

ShangHai International Holding Corp. GmbH (Europe)

Address: Eiffestrasse 80, 20537 Hamburg Germany

**Shenzhen Kingyield Technology Co.,Ltd.**

Section C1-2, FuHai Industrial Zone

Fuhai Road FuYong Town

Bao'an District

518103 ShenZhen City, GuangDong

PEOPLE'S REPUBLIC OF CHINA

## INDEX

---

INTRODUCTION.....	19
SYMBOLS IN THIS MANUAL AND PRODUCT.....	19
USING AND SAFETY INSTRUCTIONS.....	20
KNOW YOUR UNIT.....	21
BATTERY INSTALLATION/REPLACEMENT.....	22
ADJUSTING SETTINGS.....	22
BEFORE TAKING A MESSUREMENT.....	24
PROPER WAY OF MEASUREMENT.....	24
TAKE A MEASUREMENT.....	25
USE THE MEMORY FUNCTION.....	26
IIRREGULAR HEARTBEAT SYMBOL.....	27
ABOUT BLOOD PRESSURE.....	27
ERROR INDICATORS.....	29
ELIMINATE FAULT.....	30
CLEANING AND MAINTENANCE.....	30
SPECIFICATIONS.....	31
INDICATIONS FOR USE STATEMENT.....	32
GUIDANCE AND MANUFACTURE'S DECLARATION.....	32

**INTRODUCTION**

- Thank you for having chosen a blood pressure monitor ORBEGOZO.
- The unit uses the oscillometric method of blood pressure measurement. It means the unit detects the movement of your blood through your brachial artery, and converts your blood pressure into a digital reading. The unit is simple to use because a stethoscope is not needed while using an oscillometric monitor.
- The unit stores automatically 60 sets of measurement values for each user (2x60). You can read the stored data conveniently by pressing the memory button.
- The unit comes with the following components: main unit, arm cuff and instruction manual

**SYMBOLS IN THIS MANUAL AND PRODUCT**

	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.		
	Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the equipment or other property.		
	Type B applied part		Transport package should be kept away from rain.
	Authorised representative in the European community		
	Refer to instruction manual/booklet		
	Manufacturer		
	Serial number		
	The marking of electrical and electronics devices according to Directive 2002/96/EC. The device, accessories and the packaging have to be disposed of waste correctly at the end of the usage. Please follow Local Ordinances or Regulations for disposal.		



## USING AND SAFETY INSTRUCTIONS


To assure the correct use of the product basic safety measures should always be followed including the precautions listed below:

- Read all information in the instruction manual and any other literature included in the box before using the unit.
- Contact your physician for specific information about your blood pressure. Self-diagnosis and treatment using measured results may be dangerous. Follow the instructions of your healthcare provider.
- Operate the unit only as intended. Don't use for any other purpose.
- The unit is intended to use in measuring blood pressure and pulse rate for adult only, not recommended for neonatal baby at home or medical centre.
- Do not disassemble or attempt to repair the unit or components.
- Do not use the equipment where flammable gas (such as anesthetic gas, oxygen or hydrogen) or flammable liquid (such as alcohol) are present.
- Do not use a cellular phone near the unit. It may result in operational failure.
- Indicating that too frequent measurements can cause injury to the patient due to blood flow interference.
- Remove the batteries if the unit will not be used for three months or more.
- The application of the cuff over a wound can cause further injury.
- Warning notices:
  - Regarding the effect of blood flow interference and resulting harmful injury to the patient caused by continuous cuff pressure due to connection tubing kinking
  - Regarding the application of the cuff and its pressurization on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt is present because of temporary interference to blood flow and could result in injury to the patient.
  - Regarding the application of the cuff and its pressurization on the arm on the side of a mastectomy.
  - Regarding the information that pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring ME Equipment on same limb
  - Regarding the need to check (for example, by observation of the limb concerned) that operation of the automated sphygmomanometer does not result in prolonged impairment of the circulation of the blood of the patient.


KNOW YOUR UNIT

- 2 AAA (LR03) Alkaline Batteries


- Storage Case

## BATTERY INSTALLATION/REPLACEMENT

1. Remove the battery cover from the battery compartment.
2. Install 2 AAA size batteries so the + (positive) and - (negative) polarities match the polarities of the battery compartment as indicated.
3. Replace the battery cover.

## ATENCIÓN

- Do not use rechargeable batteries.
- Use two identical batteries type AAA.
- Replace both of the batteries if the symbol of low battery appears on the display.
- Never leave any low battery in the compartment since they may leak and cause damage to the unit.
- Disposal of used batteries should be carried out in accordance with the national regulations for the disposal of batteries.

## ADJUSTING SETTINGS

## 1. To set year

All information will appear when the batteries are installed at the first time. After 2 seconds the year flashes.


Press the "MEM" button to increase one year.

Press the "SET" button to decrease one year.

Long press the button and the value quickly changes.


The year can be set between 2011 and 2030, and it will return to 2011 when the display of year up to 2030.

2. To set date and time


Press the "MEM" button to change the month and press the "SET" button and continue with next step.

Then set the day, the hour and the minute by using the same way. Continue with step 3.


3. To set unit

- Then press the "MEM" button again to change the current unit.
- Press the " START-STOP" button to save the current setting and turn off.

4. To set user

With the device turned off, the mark of a user will appear when pressing the SET button. Then press SET again to change the current user. Finally, press START-STOP button to confirm user and turn off.



**BEFORE TAKING A MEASUREMENT**

Follow these instructions before taking a measurement:

1. Avoid eating, drinking alcohol, smoking, exercising and bathing for 30 minutes before taking a measurement.
2. Stress raises blood pressure. Avoid taking measurements during stressful times.
3. The cuff can be applied to your left or right wrist. Always measure your blood pressure on the same Arm.
4. Measurements should be taken in a quiet place.
5. Position the unit at heart level throughout the measurement.
6. Remain still and do not talk during the measurement.
7. Keep a record of your blood pressure and pulse readings for your physician. A single measurement does not provide an accurate indication of your true blood pressure. You need to take and record several readings over a period of time. Try to measure your blood pressure at the same time each day for consistency.
8. Wait 30-60 seconds between measurements. The wait time allows the arteries to return to the condition prior to taking the blood pressure measurement. You may need to increase the wait time depending on your individual physiological characteristics.


**PROPER WAY OF MEASUREMENT**

<p>Remove all clothes from your arm allowing the cuff to fit directly on the skin.</p>	
<p>Wrap the cuff tightly around your wrist by using the Velcro strip.</p>	
<p>Sit in a chair with your feet flat on the floor.</p>	


**TAKE A MEASUREMENT**

1. Press the START-STOP button. All display symbols appear on the screen. The cuff starts to inflate automatically.


2. Inflation stop automatically and the measurement is started. The heartbeat symbol flashes once a pulse is

detected.


3. The blood pressure and pulse rate are displayed when the measurement completed. The cuff is deflated automatically, and all of the measurement results are stored in the memory with date and time. The "IHB" symbol will be displayed if irregular heartbeat is detected.

4. Press the "START-STOP" button to turn off the unit. The unit will automatically turn off after 2 minutes if you have forgotten to turn it off.

NOTE: To stop the inflation or measurement, press the "START-STOP" button in anytime. The unit will stop inflating and then release and turn off.


The unit stores the blood pressure and pulse rate in the memory after completing a measurement every time, up to 2x60 sets. The earliest record is deleted to set the latest measurement taken.

USE THE MEMORY FUNCTION


1. To read the average value

Press the "MEM" button. The average value symbol will display.

2. To read the measurement value

Press the "MEM" button to view the latest measurement value with the measuring date and time.  
 Press the "SET" button to view the oldest measurement value with the measuring date and time.

The 1st latest measurement value

The 2nd latest measurement value

The 2x60th latest measurement value

3. To delete the memory value

In the memory reading mode, press the "MEM" button for 4-6 seconds the unit will display à **dEL no**

Press the "MEM" button again and LCD will display à **dEL YES**


If you are sure all the data of corresponding user can be deleted, press the "START-STOP" button. All the corresponding user's data will be deleted and the unit will turn off.

Note: if you have set **DEL YES**, but do not press a button for 2 minutes, the product will shut off automatically, but the corresponding user's blood pressure data will not be deleted. The delete operation must be confirmed by user.

## IRREGULAR HEARTBEAT SYMBOL

When the unit detects an irregular rhythm two or more times during the measurement, the irregular heartbeat Symbol will appear on the display with the measurement values.


An irregular heartbeat rhythm is defined as a rhythm that is more than 25% slower or 25% faster from the average rhythm detected while the monitor is measuring the systolic blood pressure and the diastolic blood pressure. If the Irregular Heartbeat Symbol (IHB) displays with your measurement results, we recommend you consult your physician. Follow the directions of your doctor.

## ABOUT BLOOD PRESSURE

### Blood Circulation


The blood circulation is responsible for supplying the body with oxygen. Blood pressure is the pressure exerted on the arteries. The systolic blood pressure value (higher pressure or top value) represents the blood pressure produced by contraction of the heart muscle. The diastolic blood pressure value (lower pressure or lower value) represents the blood pressure produced by relaxation of the heart muscle.

### Health and blood pressure

The incidence of hypertension increases with age. In addition, a lack of exercise, excess body fat and high levels of cholesterol (LDL), which sticks to the inside in blood vessels, reduces elasticity of these vessels. Hypertension accelerates arterial sclerosis which can lead to very serious conditions such as stroke and myocardial infarction. For these reasons it is very important to know whether the blood pressure is within a healthy range.


Blood pressure fluctuates from minute to minute, throughout the day. Therefore it is essential to take regular measurements to help you identify an average blood pressure.

**Classification of blood pressure**

After each measurement is completed, the LCD display will show your position automatically on the six segments of the bar indicator which corresponds to World Health Organization (WHO) Blood Pressure Indicator.

\*Note! When a person’s systolic and diastolic pressures fall into different categories, the higher category should apply.

### Symptoms of high blood pressure

High blood pressure can go unnoticed for a long time, since it does not cause noticeable symptoms. The following are all possible causes of abnormally high blood pressure:

- Overweight
- High cholesterol level
- Smoking
- Excessive alcohol consumption
- Stress and emotional upset
- Excessive consumption of salt
- Lack of physical exercise
- Genetic/hereditary predisposition
- Underlying illnesses, such as kidney disorders or endocrine disturbance

### Treatment of high blood pressure

If your blood pressure reaches upper values of 140-160 mmHg and lower values of 90-95 mmHg in repeated measurements over several days, you should consult your doctor for detailed medical examination. You can assist the treatment prescribed by your doctor in the following ways:

- Lose weight and lower you cholesterol level
- Reduce the consumption of alcohol and salt
- Stop smoking
- Take regular exercise
- Monitor your blood pressure

### ERROR INDICATORS


SYMBOL	CAUSE	CORRECTION
	The course of inflating appears error.	Wrap the cuff correctly and tightly.
		Ensure the connection is close between the air plug and the unit.
		Inflate over again after ensuring.
	When measurement fails.	Do not move your arm and body, and keep quiet.
Measure over again according to correct way.		
	When the batteries power are too low.	Replace all of the worn batteries with new ones.

**ELIMINATE FAULT**

PROBLEM	CAUSES AND SOLUTIONS
No power.	Replace all the worn batteries with 4 new ones.
No display appears on the display screen.	Check the battery installation for proper placement of the batteries polarities.
Measurement values appear too high or too low.	Blood pressure varies constantly. Many factors including stress, time of day, how you wrap the cuff, may affect your blood pressure. Review the sections "Proper Use of the Arm Cuff" and "Take a Measurement".

**CLEANING AND MAINTENANCE**

- Do not disassemble the unit
  - Do not fold the cuff or immerse it in water
  - Do not use volatile liquids to clean the main unit
  - Do not use petrol, thinners or similar solvents to clean the cuff
  - Do not subject the unit to strong shocks or vibrations
  - Do not subject the main unit and the cuff to extreme temperatures, humidity, moisture or direct sunlight.
  - Do not carry out repairs of any kind yourself. If a defect occurs, please contact your distributor.
- à If needed, please use ultraviolet radiation of disinfection.
- The unit should be cleaned with a soft, dry cloth.
  - Use a soft, moistened cloth and soap to clean the arm cuff
  - Substitution of a cuff different from that supplied might result in measurement error.

The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life. We recommend that the performance be checked every 2 years and after maintenance and repair, by utilizing the monometer mode and verifying the accuracy of the manometer at least at 50 mmHg and 200 mmHg.

## SPECIFICATIONS

Model No.		TES 3650
Display		LCD Digital Display
Measuring Principle		Oscillometric Method
Measurable Circumference of Wrist		13.5~21.5cm (5.3"~8.5")
Electric Shock Protection		Internal power supply appliance type BF
Measurement Range	Pressure	0 to 299 mmHg (0~39.9 kPa)
	Pulse	40 to 180 Pulses/min
Accuracy	Pressure	$\pm 3$ mmHg (0.4kPa)
	Pulse	$\pm 5\%$ of reading
LCD Indication	Pressure	3 Digits Display of mmHg or kPa
	Pulse	3 Digits Display
	Symbol	Memory/IHB/Heartbeat/Low Battery, etc.

Inflation	Automatic by Internal Pump
Rapid Air Release	Automatic by Air Valve
Memory Function	Memory 2x60 Sets Measurement Values
Power Source	2 AAA (LR03) Alkaline Batteries
Automatic Turn Off	Within 2 Minutes
Gross Weight	Approx. 90g (Batteries not included)
Device Life	2 Years

Operating Environment	Temperature	10~40 °C (50 ~ 104 °F)
	Humidity	15 ~ 90% RH (Noncondensing)
	Barometric Pressure	80~105 kPa
Transport Storage Environment	Temperature	-20~60 °C (-4 ~ 140 °F)
	Humidity	10 ~ 95% RH (Noncondensing)

Note: Subject to technical modification without prior notice.

## INDICATIONS FOR USE STATEMENT

- The unit is intended to be used by adults to measure blood pressure and pulse rate.
- The unit is intended of using in only adult population, not applied to the other populations such as neonatal baby.
- It cannot be used while the wrist (arm) has bleeding or wound to avoid the blood flowing from the wound in pressuring.
- The device, part and batteries have to be disposed of waste correctly at the end of the usage.
- Please follow local ordinances or regulations for disposal.
- Applied part: cuff
- Protection class: internally powered equipment.
- Applied part type: type BF.
- Moisture protection: IPX0, continue operation.
  - Altitude <2000 m.
  - Overtoltage: II
  - Pollution degree: 2
- The risk of patient and user can be lowered to acceptable level
- The unit might not meet its performance specifications if stored or used outside the following specified temperature and humidity ranges.
- The unit satisfies the requirements of IEC60601-1 Medical electrical equipment; IEC60601-1-2:2007: Electromagnetic compatibility-Requirements and tests; EN1060-1; Non-invasive sphygmomanometers-Part 1: requirements and test methods for non-automated measurement type; EN1060-3: Non-invasive sphygmomanometers – Part 3: supplementary requirements for electro-mechanical blood pressure measuring systems.

## GUIDANCE AND MANUFACTURE'S DECLARATION

ELECTROMAGNETIC EMISSION		
The Digital Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer of the user of the Digital Blood Pressure Monitor should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
Radio frequency (RF) emissions CISPR 11	Group 1	The Digital Blood Pressure Monitor use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Digital Blood Pressure Monitor is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

ELECTROMAGNETIC IMMUNITY			
The Digital Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of Digital Blood Pressure Monitor should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 V/m	3 V/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz a 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Digital Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d=1.167\sqrt{P}$ <p>80 MHz to 800 MHz</p> $d=2.333\sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Digital Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the Digital Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Digital Blood Pressure Monitor.			
<sup>b</sup> Over the frequency range 150 KHz to 80MHz, field strengths should be less than 3V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the Digital Blood Pressure Monitor.

Pressure Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Digital Blood Pressure Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Digital Blood Pressure Monitor as recommended below, according to the maximum output power of the communications equipment.

Output power of transmitter	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d=1.167\sqrt{P}$	80 MHz to 800 MHz $d=2.333\sqrt{P}$
0.01	0.117	0.233
0.1	0.369	0.738
1	1.167	2.333
10	3.689	7.379
100	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Disposal of old electrical appliances.

The European directive 2002/96/EC on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment. The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected. Consumers should contact their local authority or retailer for information concerning the correct disposal of their old appliance.

DECLARATION OF CONFORMITY: This device complies with the requirements of the Low Voltage Directive 2006/95/EC and the requirements of the EMC directive 2004/108/EC.

European representative name:  
ShangHai International Holding Corp. GmbH (Europe)  
Address: Eiffestrasse 80, 20537 Hamburg Germany

Shenzhen Kingyield Technology Co.,Ltd.  
Section C1-2, FuHai Industrial Zone  
Fuhai Road FuYong Town  
Bao'an District  
518103 ShenZhen City, GuangDong  
PEOPLE'S REPUBLIC OF CHINA

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